isogait

Hydrotherapy Treadmill

Underwater treadmill for physiotherapy and rehabilitation




Designed for supervised walking and movement training, with adjustable treadmill speed, automatic water-level control and programmable session presets.

7-inch	0.05–10 km/h	Up to 37°C
Japanese touchscreen	Treadmill speed	Temperature control

isogait

Simple, precise control



01

Japanese touchscreen

7-inch operator interface for session settings.

02

Glass water chamber

A clear view of movement inside the chamber.

03

Internal handrails

A handhold during supervised movement.

Adjustable treadmill speed

Set walking speed from 0.05 to 10 km/h.

Water level and temperature

Automatic water-level control; temperature control up to 37°C.

Programmable presets

Save session settings for repeat use.

isogait

Built for water-based movement training

For supervised use in physiotherapy clinics, rehabilitation facilities and sports recovery settings.

Movement in water

Buoyancy helps support movement; water provides resistance.

Specification	Detail
Construction	SS 304 stainless steel and glass
Water chamber dimensions	W 32 × L 54 × H 53 in
Overall machine dimensions	W 33 × L 60 × H 66 in
Machine sizing	Customisable
Treadmill speed	0.05–10 km/h
Touchscreen	7-inch Japanese touchscreen
Controller	Japanese controller
Water temperature	Automatic control up to 37°C
Water level	Automatic control
Filtration	Sand and UV filtration
Session settings	Programmable presets
Bearings	Corrosion-resistant
Optional features	Jacuzzi and lighting

isogait

Plan your hydrotherapy setup

Contact isogait for product details, configuration options and a quotation.



Enquire with isogait

Hydrotherapy treadmill manufacturer

Call / WhatsApp

+91 77386 09303

Email

sales@isogait.com

Website

www.isogait.com